

**Product Information** 

ISO Class 3 - 5

Cleanroom Class 1 – 100 EU Grade A/B/C/D

# Sterile SatPax® HA

Sterile Pre-wetted Sealed Edge Cleanroom Laundered 100% Polyester Knit Wiper

Sterile SatPax® HA combines ValuSeal® HA ultrasonically sealed edge cleanroom laundered wipers composed of 100% continuous filament polyester knit fabric with a high saturation level of 70% IPA and 30% DI water. This pre-wetted format provides a cost-effective and easy-to-use solution versus traditional bulk handling of solvents, maintenance of squirt bottles and inconsistent wetting and cleaning associated with wetting a dry wiper.



#### Other Sterile Pre-wetted wipers

- Sterile SatPax® 670
- Sterile SatPax® 1000
- Sterile SatPax® 1200

## **Key Attributes**

- 100% continuous filament polyester knit
- Ultrasonically sealed edge for reduced fibre contamination
- · Laundered and packaged in Berkshire's ISO Class 4 cleanroom
- Pre-wetted with consistent 70% IPA/ 30% DI Water to a high saturation level
- Re-sealable solvent resistant packaging
- · Wipers in j-folded configuration for single withdrawal
- Gamma irradiated and sterile validated to a 10<sup>-6</sup> Sterility Assurance Level per AAMI guidelines
- Certificate of Sterility available to down load on the web by lot number **24/7**. Details: Expiration Date, Radiation Dosage, Lot Information, Meets cGMP requirements for traceability

#### **Benefits**

- Critically low fibres, particles, ions and extractables combined with high abrasion resistance
- Passes LAL test for endotoxins <20 EU/device</li>
- Reduces alcohol usage and preparation / handling costs
- Reduces VOC emissions
- Increases cleaning efficiency
- Increases cleaning protocol consistency

### **Environmental**

Registered under REACH

#### **Applications**

- Designed for use in ISO Class 3 and higher sterile cleanroom environments
- Designed for use in wet cleaning of critical surfaces where control of flammable solvents and flammable solvent
  concentrations is required
- Final cleaning of surfaces or products prior to manufacturing or packaging
- High saturation level is ideal for removing cleaning and disinfecting residues in sterile environments

#### **Saturation Levels**

The amount of solution contained in each wiper can be varied according to customer requirements. Higher saturation levels apply more solution to the surface during cleaning.

## Validated Sterile

ANSI/AAMI/ISO procedures have been followed for verification of irradiation dose with quarterly audits being performed to monitor bioburden levels and dose. The Sterility Assurance Level (SAL) or the probability of a viable organism being present on a product unit after sterilization is a SAL of  $10^{-6}$  (which is the chance of 1 non-sterile unit in 1,000,000).

## www.berkshire.uk.com

Contact: Tel + 44 1953 562800 enquiries@berkshire.uk.com

America	Tel 1 413 528 2602	info@berkshire.com				
Europe	Tel + 44 1953 562800	enquiries@berkshire.uk.com				
SE Asia	Tel 65 6252 4313	enquiries@berkshire.com.sg				
Japan	Tel 81 3 4530 9883	master@berkshire.co.jp				



## **Technical Data (In Dry State)**

Attribute		Units	Value	Test Method		
Basis Weight		g/m²	136	TAPPI T-410		
Caliper	_	μm	435	TAPPI T-411		
Fibres	≥100µm	fibres/cm²	0.080	IEST-RP-CC004.4 Sec 7.1.3/Sec 7.2.2 modified		
Particles	<u>&gt;</u> 0.5μm	x10³/cm²	0.31	IEST-RP-CC004.4 Sec 7.1.3/Sec 7.2.1 modified		
Sorbency	Capacity	mL/m²	567	IEST-RP-CC004.4 Sec 9.1 / Sec 9.2 modified		
	Efficiency	mL/g	4.2			
	Rate	seconds	1			
Non-Volatile Residue	DI Water	g/m²	0.0042	IEST-RP-CC004.4 Sec 8.1.2		
	IPA	g/m²	0.0092			
Ions	Na <sup>+</sup>	ppm	0.41	IEST-RP-CC004.4 Sec 8.2.2		
	K <sup>+</sup>	ppm	0.048			
	Ca <sup>++</sup>	ppm	0.035			
	Mg <sup>++</sup>	ppm	0.0085			
	CI <sup>-</sup>	ppm	0.21			

#### Notes:

- Technical data represented in this table are typical values at the time of publication. These should not be used as
  product specifications.
- Due to differences in test methods applied and equipment utilised by different wiper manufacturers, valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions
- Third party testing can be performed upon request

## **Order Information:**

Product	Number	Size	Shts/pk	Pks/cs	IPA/DI	Saturation	VOC % by Weight	Style
Sterile SatPax® HA	SSPHA00112	9x9" (23x23cm)	30	12	70/30	60%	42%	J-fold

## Other Berkshire products



Wipers



**Glove Liners** 



**Mop Systems** 



**Documentation Systems** 



**Face Masks** 



Swabs